Gebrauchsanleitung **Armsysteme**  
Instructions for Use **Arm Systems**  
Istruzioni per l'uso **Il sistema per braccio**  
Manuel d'utilisation **Le système de bras**  
Instrucciones de uso **El sistema de brazo**

**HERCULES**<sup>®</sup>  
BY GEOMED

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## INTENDED USE

- The HERCULES® Arm System is a table-mounted instrument and endoscope holder system intended for the following uses in accordance with the chosen device holder:
  - When used with a Universal Holder, the HERCULES® Arm System is intended to hold retractors and other surgical instruments with a diameter of 2.7mm to 16mm during surgical procedures.
  - When used with a Universal Endoscope Holder, the HERCULES® Arm System is intended for use by surgeons to hold endoscopes and arthroscopes with a diameter of 2.7mm to 10mm during general diagnostic and therapeutic procedures. The device is also intended for use by qualified surgeons for holding endoscopes during diagnostic and therapeutic neurologic procedures.
- The Clamp Socket is used to fasten the HERCULES® Arm System to a table rail or tube with a height of 25-32 mm and a thickness of 8-12 mm.
- The HERCULES® Arm System is available in multiple sizes that differ in holding power. See table at the end of these instructions for use.

## USE OF DEVICE

- The HERCULES® Arm System may only be used for its intended purpose by properly trained and qualified personnel.
- Proper selection of the device for each application, proper surgical technique for the use of all instruments and ensuring that OR personnel is adequately trained and briefed is the responsibility of the surgeon.
- Incorrect handling and care as well as misuse or use for a task other than that for which the HERCULES® Arm System is intended could result in a damaged or broken device or in unsatisfactory performance, in addition to causing premature wear and thus reducing the useful life of the product.
- Failure to follow these instructions will render any warranty claim null and void.

## WARNINGS

- "Lubrication of the ball-and-socket joints may cause loss of holding power and collapse of the device. NEVER lubricate the ball-and-socket joints. Only the central joint with the adjustment knob may be lubricated."
- Insufficient tightening of the HERCULES® Arm System, including the Clamp Socket, may result in inadvertent movement. The device and/or clamp socket could move or drop, potentially resulting in injuries.
- When twisting the adjustment knob counterclockwise to release the joints, never turn it past the point of resistance.
- "Never exceed the maximum holding power of the device. Do not overload it. Overloading the device may cause

the HERCULES® Arm System to break and lead to injury of the patient."

- Always grasp the device with one hand prior to adjusting its position. Use the other hand to loosen the adjustment knob on the central joint in order to avoid sudden abrupt downward movement.
- The device can transmit electric current and heat to the patient, causing electric shock and/or burns. Avoid contact between any part of the device and sources of electricity or heat.
- Do not move the table with the patient while the HERCULES® Arm System is in use.
- Always remove the device after clinical use before taking the patient off the table.
- Regularly and carefully inspect the plastic parts of insulated clamp socket for cracks, porosity and damaged surfaces. Damaged plastic parts may reduce dielectric strength.
- Release all joints prior to cleaning and sterilization. Releasing the joints is necessary to ensure adequate reprocessing of the device.
- Do **not** lock the Clamp Socket (1) without a stand (2) inside!

## PRECAUTIONS

- Set down and store products in accordance with their mechanical sensitivity to avoid damage.
- Do not combine the system components of this device with components of other manufacturers.
- Never let blood and/or other contaminants, such as body fluids, tissue or saline solution, dry on your HERCULES® Arm System. If these contaminants are allowed to dry, they may be difficult to remove and interfere with subsequent cleaning and sterilization.
- Thorough and effective cleaning is indispensable to ensure effective terminal sterilization.
- Minimize or eliminate time delays between processing steps. Delays may create conditions favorable to microbial growth or colonization and increase the challenge to each subsequent step.
- If the HERCULES® Arm System comes in contact with corrosive substances, such as blood, body fluids, saline solution or detergents, rinse the device immediately. If these substances are left to dry, they may reduce the useful life of your HERCULES® Arm System.
- For selection of an appropriate cleaning and/or disinfection agent, consult section Material Compatibility. Unsuitable products can damage product surfaces and accessories if used repeatedly in the long term.
- Never use metallic brushes or steel wool to clean the device components.
- Never oil plastic components.
- Never give a device a complete coating of oil.
- Never subject an HERCULES® Arm System to temperatures above 285°F (141°C).

## NEW INSTRUMENTS

- Upon receipt, remove all protective packaging and dispose of it in an environmentally friendly way before preparing the product for processing and clinical use. The product is supplied non-sterile.
- Thoroughly inspect your new HERCULES® Arm System for damage that may have occurred during transit. Immediately return damaged devices to the manufacturer

## DEVICE ASSEMBLY FOR USE

### Anchoring the HERCULES® Arm System

- Secure the Clamp Socket (1) to the operating room table side rail or tube by turning the rocker arm counterclockwise
- Insert the stand (2) of the HERCULES® Arm System into the perforation of the Clamp Socket.
- Turn the lever (3) of the Clamp Socket clockwise to fix the HERCULES® Arm System at the desired height.
- Do **not** lock the Clamp Socket (1) without a stand (2) inside!

- Grasp the distal segment of the HERCULES® Arm System firmly with one hand.
- Use your other hand to turn the adjustment knob on the central joint (5) counterclockwise, releasing the three joints.
- Do not turn the knob past the point of resistance!
- The joints can now be moved freely. Bring the arm into

the desired position.

- Lock the joints of the arm by turning the adjustment knob on the central joint clockwise at **least 1½ turns**. All joints will lock simultaneously.
- You can now let go of the distal segment of the HERCULES® Arm System.

### Quick-release Adapter System

- Insert the plug located at the base of the retractor/surgical instrument or endoscope holder into the Quick-release adapter mounted on the tip of the support arm and turn the holder slightly until you hear a distinct click.
- If attached correctly, it should be impossible to turn the holder any further or to pull it out.

### Universal Holder

- Place the retractor or instrument shaft on the V-shaped notch located on the bottom inside surface of the holder
- Turn the adjustment screw clockwise tightly against the retractor or instrument shaft to secure it in place.

### Endoscope Holder

- Rotate the thumb wheel on the Universal Endoscope Holder counterclockwise to open the holder. Insert the endoscope into the desired position.
- Secure the endoscope in place by rotating the thumb wheel clockwise.

### Final Positioning of the HERCULES® Arm System

- The HERCULES® Arm System can be adjusted freely within its articulating radius. The adjustment knob on the central joint is designed to simultaneously lock all the joints of the system into any desired position.
- The device may vibrate slightly when touched. This is normal due to its length and design.

### DEVICE DISASSEMBLY AFTER USE

- Grasp the instrument or endoscope with one hand and turn the adjustment screw on the holder counterclockwise in order to release and extract the instrument or endoscope from the holder.
- If you are using the HERCULES® universal holder, use one hand to depress the knob located on the Quick Release Adapter mounted on the tip of the support arm and extract the retractor or endoscope holder by pulling with your other hand.
- Firmly hold the stand of the HERCULES® Arm System with one hand. Turn the lever on the Clamp Socket counterclockwise to release the HERCULES® Arm System. Vertically lift it out of the Clamp Socket and lay it horizontally on a flat surface.
- Release the Clamp Socket from the table side rail or tube.
- The HERCULES® Arm System is now ready to be re-

processed in accordance with the instructions provided in these instructions for use.

## REPROCESSING THE HERCULES® ARM SYSTEM

### General Information

- HERCULES® Arm Systems are supplied non-sterile and must be cleaned and terminally sterilized before first use and after each reuse.
- Thorough and effective cleaning is indispensable to ensure effective sterilization.
- The cleaning and sterilization procedure provided below has been validated by the manufacturer as being capable of preparing the HERCULES® Arm System for reuse and achieving sterilization to a sterility assurance level (SAL) of 10<sup>-6</sup>.
- "To ensure safe and effective reuse of the device, it remains the responsibility of the user to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired results. This requires verification and/or validation and routine maintenance and monitoring of the process."
- All processing steps must be performed by adequately trained personnel using the appropriate equipment and materials under strict observance of the validated parameters.
- Observe nationally applicable guidelines and legal provisions as well as the hygienic regulations of your facility, particularly regarding effective inactivation of prions. Any deviation from the processing procedure described in these instructions for use must be validated by the user.
- These reprocessing instructions conform to the requirements of EN ISO 17664-1 – Processing of health care products – Information to be provided by the medical device manufacturer for the processing of medical devices.

### PRECAUTIONS

- Never let blood and/or other contaminants, such as body fluids, tissue or saline solution, dry on your HERCULES® Arm System. If these contaminants are allowed to dry, they may be difficult to remove and interfere with subsequent cleaning and sterilization.
- Thorough and effective cleaning is indispensable to ensure effective terminal sterilization.
- Minimize or eliminate time delays between processing steps. Delays may create conditions favorable to microbial growth or colonization and increase the challenge to each subsequent step.
- Follow the instructions and warnings issued by the suppliers of your cleaning and sterilization equipment.
- If the HERCULES® Arm System comes in contact with corrosive substances, such as blood, body fluids, saline solution or detergents, rinse the device immedi-

ately. If these substances are left to dry, they may reduce the useful life of your HERCULES® Arm System.

- For selection of an appropriate cleaning and/or disinfection agent, consult section Material Compatibility. Unsuitable products can damage product surfaces and accessories if used repeatedly in the long term.
- Do not expose the device to temperatures above 285°F (141°C).

## REPROCESSING INSTRUCTIONS

### During Use

• If the device is not to be used for the remainder of a procedure, prevent the drying of blood, soil and debris on its surface by wiping it clean using a disposable wipe or by placing it to soak in a basin containing purified cold water.

### Point of Use

- Gross soil and blood must be removed from the HERCULES® Arm System as soon as possible after use. To facilitate the process, lock the joints by turning the adjustment knob on the central joint clockwise ½ to ¾ turn. The device will remain in a rigid position, making it easier to handle.
  - Open the Clamp Socket completely.
  - Wipe gross soil and blood from the HERCULES® Arm System using a disposable wipe
  - Remove any visible contamination from the device using a soft brush or a clean, lint-free cloth while rinsing the device thoroughly with one to two liters running deionized water at 35°C-40°C for one to two minutes.
  - While rinsing the socket, shift the movable bracket back and forth several times.
  - Place the HERCULES® Arm System in an instrument tray or container.
  - Keep the HERCULES® Arm System wet by covering with a towel moistened with purified water or by using a foam, spray or gel product specifically intended for this use.
  - Keep the outside surface of the container clean.
  - Transport the tray to the decontamination area in an impermeable plastic bag or a container with a tight lid.
  - The HERCULES® Arm System may additionally be disinfected at point of use.
- ### Optional Disinfection at Point of Use:
- For protection of hospital personnel, the HERCULES® Arm System may be disinfected by immersion in a freshly prepared disinfectant solution at point of use according to the disinfectant manufacturer's instructions for use.
  - The disinfection procedure is intended for protection of hospital personnel only. Disinfection cannot replace the final sterilization process.
  - If you choose to disinfect the HERCULES® Arm Sys-

tem, the disinfectant used must be alcohol and aldehyde-free, compatible with instruments manufactured from stainless steel and thermoplastic and legally marketed for use with medical devices. (See "Material Compatibility.")

- Change the soaking solution once a day or more often if heavily contaminated.

### CLEANING

• Process the device as soon as is reasonably possible after use. It is recommended not to delay the main cleaning procedure for more than 2 hours.

### Pre-cleaning

- Wiping off the coarse soiling with disposable cloths.
- Brushing of the objects with soft brushes under cold running deionized water for 1 minute, while rinsing move the moving parts back and forth 3 times each.

### Cleaning Process

- The automatic cleaning instructions described in this section have been validated by the manufacturer using the mildly enzymatic alkaline cleaning agent neodisher® MediClean Forte by Dr. Weigert which has a pH range of 10.4-10.8 (2-10 ml/l, determined in deionized water at 20°C).
- If you will be using a different cleaning agent, make sure that the product has equivalent specifications and is suitable for cleaning of surgical instruments made of stainless steel and thermoplastic. A low sudsing/foaming detergent should be preferred.
- Always follow the detergent manufacturer's instructions regarding the concentration of the solution, water quality and temperature as well as soaking and rinsing times.

### Automated Cleaning and Disinfection

- Release all joints of the arm component to ensure adequate cleaning and disinfection of inner joint surfaces.
- Place the HERCULES® Arm System components in an appropriate instrument tray or basket, making sure that the devices do not touch.
- Start the cleaning program. The following minimum parameters are recommended.
  - 1 minute prewash with cold tap water
  - 3 minutes precleaning with cold tap water
  - 5-10 minutes main cleaning with a 0.5% cleaning agent and deionized water at 131°F (55°C)
  - 2 minutes final rinse in cold, deionized water
  - 5 minutes thermal disinfection with deionized water at 194°F (90°C)
- It is recommended to dry the components of the HERCULES® Arm System hot air dry with filtered air for 15-30 minutes at 230°F (110°C) to dry. Drying time is shown as a range because it depends on the

load size placed in the washer/disinfectant.

- After completion of the program, remove the devices from the washer.
- If the devices are not completely dry, supplement the drying step by air-drying in a clean area or by drying with a clean absorbent, non-shedding and lint-free cloth.
- Clean, filtered compressed air may be used to remove moisture from difficult-to-access areas.
- All surfaces must be completely dry before sterilization.
- Inspect the HERCULES® Arm System components for damage.

## INSPECTION AND FUNCTION CHECK

- After cleaning, perform the following inspection and function check.
- If damaged components are identified, place them aside, making sure to mark them as non-conforming.
- The symbol ☒ denotes an action that must be taken by the user.
- The symbol ☑ denotes the reaction produced by such an action.
  - Example:
    - ☒ Release the joint.
    - ☑ The joint has been released.

### No. 1: Device dry?

Action:

- ☑ Dry all openings with compressed air as needed.
- ☑ Device dry. Continue to step #2.

### No. 2: Device clean?

Action:

- ☑ Visually inspect the device for cleanliness.
- ☑ No residues remain. Continue to step #3.

☑ Residue is visible on the device.

- ☑ Repeat the cleaning procedure and restart the inspection and function check from step #1.

### No. 3: Device undamaged?

- ☑ Visually inspect the device for corrosion, for damaged, splintered, chipped, cracked or rough surfaces and for sharp edges.
- ☑ No damage found. Continue to step #4.

☑ Damage found.

- ☑ Send product to manufacturer for repair. See "Repair" and continue to step #6.

### No. 4: Rods and joints

- ☑ Rods are straight and firmly connected to central joint and ball-and-socket joints.
- ☑ Continue to step #5.

- ☒ Rods are bent or not firmly connected to central joint and/or ball-and-socket joints.
- ☒ Send product to manufacturer for repair. See "Repair" and continue to step #6.

#### No. 5: Joint locking and release mechanism

- ☒ Lock joints by rotating the adjustment knob located at the central joint clockwise.
- ☒ All 3 (three) joints lock correctly.
- ☒ If joints do not lock correctly, stop and send product to manufacturer for repair. See "Repair" and continue to step #6.
  
- ☒ Release the correctly locked joints by rotating the adjustment knob located at the central joint counterclockwise.
- ☒ All three joints are released and move easily. Continue to step #6.
  
- ☒ Joints are released, but do not move easily.
- ☒ Lubricate the central joint with a small amount of oil in accordance with the instructions provided in section "Lubrication".  
**WARNING: Only the central joint with the adjustment knob may be lubricated. NEVER lubricate the ball-and-socket joints.**
- ☒ All three joints are released and move easily. Continue to step #6.
  
- ☒ Despite lubrication of the central joint, the joints do not move easily.
- ☒ Send product to the manufacturer for repair. See "Repair" and continue to step #6.

#### No. 6: Clamp Socket

- ☒ Carefully inspect the plastic parts of the insulated clamp socket for cracks, porosity and damaged surfaces. Damaged plastic parts may reduce dielectric strength.
- ☒ No damage found. Continue to step #7.
  
- ☒ Plastic parts of clamp socket are cracked, porous or damaged.
- ☒ Send Clamp Socket to manufacturer for repair. See "Repair."

#### No. 7: End of inspection procedure

- Immediately proceed to device lubrication, if applicable, and packaging for sterilization.

#### LUBRICATION

#### WARNINGS

- Lubrication of the two **ball-and-socket** joints may cause loss of holding power and collapse of the

- device. NEVER lubricate the ball-and-socket joints
- Only the **central joint** with the adjustment knob may be lubricated.

#### PRECAUTIONS

- Lubrication is not necessary on a routine basis, but moving parts made of stainless steel may be lubricated if needed.
- For lubrication, only use instrument oil of proven biocompatibility suitable for steam sterilization at 270°F (132°C).

#### Instructions

- Apply a small amount of oil to the joint.
- Move the joint to distribute the oil.

#### STERILIZATION

##### Packaging

- Release all joints of the arm component to ensure adequate sterilization of inner joint surfaces.
- Double-wrap the device components individually in single-use sterilization pouches that are approved for steam sterilization according to EN ISO 11607.
- If you are located in the United States, make sure to use packaging that has been FDA-cleared for pre-vacuum steam sterilization at 270°F with a holding time of 4 minutes, 20 minutes drying time and double packaging.

##### Sterilization

- Use a validated, properly maintained and calibrated steam sterilizer. Follow the sterilizer manufacturer instructions with regard to load configuration and equipment operation. Make sure that no rust, dirt or other foreign or corrosive substances are present in the sterilization system, which could contaminate the steam and corrode instruments, leave surface deposits or cause other damage.
- Place the wrapped HERCULES® Arm System components in the sterilizer.
- The sterilization was validated using the following minimum process parameters:

Cycle Type	Temperature	Holding Time	Drying Time
Prevacuum	270°F (132°C)	4 minutes	20 minutes

- This validation covers all standard steam sterilization processes conducted at 134°C with a minimum holding time of 4 minutes. A validated procedure is recommended for sterilization, such as the use of a sterilizer compliant with EN 285 and validation in accordance with EN ISO 17665.

- Strictly adhere to the validated sterilization parameters recommended below, including drying time. Thorough drying of the HERCULES® Arm System after sterilization is essential to inhibit any subsequent growth of waterborne microbes.
- Remove the wrapped device from the sterilizer.
- Inspect packaging to ensure that it is dry. Items with wet packaging are considered contaminated.
- Repackage and resterilize any HERCULES® Arm System with wet packaging.

#### STORAGE

- Handling, transport and storage of all devices must be managed with care so as not to affect the devices in any way.
- Set down and store products in accordance with their mechanical sensitivity to avoid damage.
- Store new or unused devices in a clean, dry and protected place.
- Store the wrapped, sterilized HERCULES® Arm System in a designated clean and dry limited-access area that is well-ventilated and provides protection from light, UV radiation, dust, moisture, insects, vermin and extreme temperature and humidity.
- Do not store sterile products in the same space or near aggressive substances, such as alcohols, acids, bases, solvents or disinfectants.
- "Inspect every package before use to ensure that the sterile barrier (sterilization pouch) is not torn or perforated and does not show signs of moisture or appear to be tampered with. If any of these conditions are present, the contents must be considered non-sterile. In this case, repeat the entire cleaning and sterilization process according to these instructions prior to using the device."

#### ADDITIONAL INFORMATION

##### Reusability and Product Lifetime

- Processing has little effect on the HERCULES® Arm System, if instructions are followed with care and the HERCULES® Arm System is not damaged or contaminated.
- Product lifetime may vary and depends on frequency of use, care and handling. Improper, ineffective and insufficient processing and maintenance can greatly reduce the lifetime of the device.
- The HERCULES® Arm System may be used until it exhibits any one of the following end-of-life indicators:
  - Corrosion
  - Nicks, dents, chips, or cracks
  - Sharp edges or rough surfaces
- Discontinue use and dispose of the HERCULES® Arm System if any of these indicators is identified.
- End of useful instrument life is generally determined by wear or damage incurred during use.

- The arm systems are tested for 1000 reconditioning cycles. If clear signs of wear or damage are visible due to heavy use, the products must be disposed of prematurely.

#### Material Compatibility

- Do not use detergents or disinfectants containing
  - Acids (pH <5)
  - Oxidizing acids
  - Bases (pH >11)
  - Organic solvents
  - Alcoholic disinfection agents
  - Benzene, phenol or ammonia
  - Halogens, halogenated hydrocarbons or sodium chloride (at high concentrations)
  - Oxidants, peroxides or hypochlorite
  - Chlorine
  - Fluorine
  - Aldehyde.

#### Repair

- The HERCULES® Arm System cannot be repaired by third parties.
- If a HERCULES® Arm System requires repair, clean and sterilize the device according these instructions for use and return it to the manufacturer.

#### Return of Used Devices

- If a device needs to be returned to the manufacturer for any reason after clinical use, it must be cleaned and sterilized according to these instructions for use prior to shipment.

#### Disposal of Device

- Dispose of irreparably damaged HERCULES® Arm Systems in accordance with established hospital protocol for disposal of surgical waste.

#### Behavior in case of serious incidents

- All incidents that occur in connection with the product must be reported to the manufacturer and the responsible national authority.

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**Medizin-Technik GmbH & Co. KG**